



## PRACTICAL TRAINING CHECK LIST

**When applying LUCAS™2 the following activities should be performed in the following order:**

### Work in teams of two. Colleague starts with manual CPR.

**Complete    Incomplete**

1. Open the bag. Push ON/OFF for 1 second to start self-test and power up LUCAS™2.

☐
☐

2. Stop CPR and place the Back Plate.

☐
☐

### Resume manual CPR.

3. Take LUCAS™2 from the bag, extend the legs and make sure the claw locks are open - pull up once and then let go of the release rings.

☐
☐

### Clock starts.

4. Connect the Upper Part to the Back Plate – start on the side closest to you. Listen for click. Pull up once to ensure attachment.

☐
☐

5. Position the Suction Cup immediately above the end of the sternum and make sure it is centered over the chest. Adjust if necessary.

☐
☐

6. Push the Suction Cup down with two fingers (make sure it is in the ADJUST mode). The pressure pad inside Suction Cup should touch patient's chest. Push PAUSE to lock start position.

☐
☐

7. Push ACTIVE (continuous) or ACTIVE (30:2) to start LUCAS™2 compressions.

☐
☐

### Clock stops.

8. Attach the Stabilization Strap and fasten the manikins arms to the patient straps, if possible.

☐
☐

Ready!

**All steps should successfully be completed and interruptions of manual CPR should be less than 20 sec.**

☐
☐

**TRAIN REGULARLY TO KEEP UP PROFICIENCY!**

Interruption of manual CPR. Max 20 sec.